


Exploring Therapeutic Potential: Insights from Drug Penetration and Targeting using Nilogen's Ex-vivo 3D-EXpress Tumoroid Platform

Nilogen Oncosystems Tampa FL 33612



Background




Figure 1: Drug Discovery Process. The schematic illustrates the sequential stages of the drug discovery process. It begins with the identification of a therapeutic target, followed by the screening of potential compounds and lead identification. Subsequent steps include lead optimization and preclinical testing to assess the efficacy and safety of candidate molecules. Promising candidates advance to clinical trials, which are conducted in phases (I, II, and III) to evaluate safety, dosage, and therapeutic efficacy in human subjects. The final stage involves regulatory review and approval before the drug is made available for clinical use.




Figure 2: Labeling therapeutics with fluorochromes aids in understanding its mechanisms. Nilogen's 3D-EX vivo models retain the tumor microenvironment found in a patient's tumor providing an optimal platform study therapeutic efficacy. Fluorescently labeled therapeutics can be used to study binding and internalization kinetics, in addition to tracking the penetration of the drugs within the tumor.

3D-Express Platform




Figure 3: Nilogen's 3D-Express tumoroid platform capabilities. The diagram illustrates the capabilities of the Nilogen 3D Express platform, which enables comprehensive analysis of drug responses using 3D tumor models. Key features include high-throughput screening, immune cell infiltration assays, cytokine profiling, and assessment of immune checkpoint activity. The platform facilitates the evaluation of therapeutic efficacy, immune modulation, and biomarker identification in a physiologically relevant 3D microenvironment, providing valuable insights for oncology research and drug development.

Results

Assay Optimization




Fig 4: Evaluating expression of Trop2 on T98G and HCT-116 cell lines. HCT-116 and T98G cell lines were stained with fluorescently labeled Trop-2 antibody and evaluated for expression of Trop2 by flow cytometry. Flow plots showing the expression of Trop2, which is enumerated in adjacent bar graph.

Fig 5: Optimizing ADC binding and internalization in spheroid model. Govitecan was labeled with BoRad Antibody labeling kit (PE-TR) and Thermofisher Zenon phloxo kit. (A) T98G and (B) HCT-116 were cultured on FocoCellate BIOFLOAT FLEX coated wells to facilitate spheroid formation and treated with labeled Govitecan for 48hrs. Representative images of two spheroids with Govitecan binding using BoRad labeling kit (PE-TR) (C) and (D) of labeled Govitecan in both the spheroids. (E) Evaluating internalization of Govitecan using Thermofisher Zenon phloxo kit in HCT-116 spheroid.

Assay Validation




Fig 6: Optimizing ADC binding and internalization in tumoroids. Endometrial tumoroid tissue were treated with therapeutics labelled to study binding (BoRad) and internalization (BoRad). Representative images showing the bound and internalized therapeutic. Graphical representation showing the cell count in each zone of penetration (μm) within the tumoroids.

How labelled therapeutics enhance research




Fig 7: Evaluating therapeutic binding to the tumoroids. (A) This figure illustrates the process of therapeutic monitoring using labelled therapeutics. Labelled therapeutic agents are introduced to the system, allowing for real-time tracking and measurement of therapeutic efficacy. The labelled molecules bind specifically to target cells or pathways, enabling visualization and quantification. (B) Representative images showing labelled therapeutic binding to RCC tumoroids. (C) Bar graph showing therapeutic bound to the cells in different RCCs and their relation to the therapeutic concentration. (D) Bar graph showing therapeutic bound to the cells in different RCCs and their relation to the therapeutic concentration and time of incubation.




Fig 8: Kinetics of therapeutic internalization in RCC tumoroids. RCC tumoroids were either left untreated or treated with two distinct labelled therapeutics to assess the rate of therapeutic internalization. The data is presented as the number of cells bound to the therapeutic at varying infiltration depths within the tumoroids, providing insight into the penetration and binding efficiency of the treatments over time.




Fig 9: Internalization of targeted vs. non-specific BTEs in RCC tumoroids. Left panel: Fluorescence microscopy images show nuclei in blue and BTEs in red. The use of labelled therapeutics allowed visualization of BTE internalization. The targeted BTE (left) shows strong internalization, while the non-specific BTE (right) shows minimal red signal, indicating weak internalization. Right panel: Bar graphs quantify the number of cells bound to the BTEs at different depths.

Summary

Evaluating Drug Mechanisms Using Labelled Therapeutics

- Efficacy**
 - Labelled therapeutics enable precise measurement of drug efficacy by tracking the therapeutic agent's activity in real time.
 - This allows researchers to assess how effectively the drug achieves its intended biological response within target tissues.
- Binding**
 - The labelled markers help in visualizing the binding process between a drug and its specific target, such as receptors or proteins.
 - This provides insights into binding affinity and helps in understanding the specificity of drug-target interactions.
- Internalization**
 - By using labelled agents, the process of drug internalization into cells can be monitored, revealing how a therapeutic enters the cellular environment.
 - This information is crucial for drugs that need to act within cells to be effective, such as targeted therapies.
- Penetration**
 - Labelled therapeutics allow the study of drug penetration through barriers like cell membranes or tissue layers.
 - This helps in evaluating the drug's ability to reach its target site within the body and contributes to optimizing formulation and delivery strategies.
- Interaction**
 - The approach provides detailed data on how drugs interact with cellular components, such as enzymes, receptors, or other molecules.
 - It helps to delineate off-target effects and elucidates the molecular pathways involved in the drug's mechanism of action.
- Mechanism of Action Evaluation**
 - Labelled therapeutics allow direct observation of drug interactions with target cells, tissues, or receptors.
 - They help in mapping the distribution and localization of drugs within the body, providing insights into pharmacokinetics and biodistribution.
 - This approach enables visualization of drug-receptor binding events, cellular uptake, and intracellular pathways.